


## Key Design Features & Instrumentation



One Tray for improved simplicity and ease of use  
Full osteotomy solution for forefoot deformity correction surgery

Foot Recon System Instrumentation		
Product Code	Description	
OS310100	FOOT RECONSTRUCTION SYSTEM STERILISATION TRAY (FRS)	1
OS310701	STAPLE HOLDER 90 DEGREE	2
OS310702	STAPLE HOLDER 26 DEGREE	3
OS310703	STAPLE IMPACTOR 90 DEGREE	4
OS310704	STAPLE IMPACTOR 26 DEGREE	5
OS322010	TWIST-OFF SCREWDRIVER	6
OS322020	META SCREW SCREWDRIVER	7
OS324000	META SCREW BONE CLAMP	8
OS325001	META SCREW BACKHAUS FORCEPS	9

3mm X 10mm - 34mm Short thread Compression screw options also available as part of this system.

### Meta Screw®

- Long and short thread options, differential pitch & self tapping
- 2 in 1 drill/counter sink size options 12mm (OS310203) 18mm (OS310204) & 22mm (OS310206) used with 1mm x 70mm trocar k-wire (OS21070H)



Instrumentation: 7 8 9

### Twist-Off

- Simple and strong fixation for Weil's osteotomy
- Reduces need to pre-drill
- Three options of twist-Off screw:
  - 1: Standard
  - 2: Wire drive
  - 3: Cancellous



Instrumentation: 6

### Varisation Staple

- Straight 90 degree 8mm & 10mm and oblique 26 degree 8mm & 10mm



Instrumentation: 2 3 4 5

# Foot Reconstruction System Product Rationale

## Product Code Reference Tables

Halgus Valgus deformity; Chevron or Scarf Osteotomy differential screws	
Product Code	Description
OS310110	3.0MM META SCREW - LONG THREAD 10MM
OS310112	3.0MM META SCREW - LONG THREAD 12MM
OS310114	3.0MM META SCREW - LONG THREAD 14MM
OS310116	3.0MM META SCREW - LONG THREAD 16MM
OS310118	3.0MM META SCREW - LONG THREAD 18MM
OS310120	3.0MM META SCREW - LONG THREAD 20MM
OS310122	3.0MM META SCREW - LONG THREAD 22MM
OS310124	3.0MM META SCREW - LONG THREAD 24MM
OS310126	3.0MM META SCREW - LONG THREAD 26MM
OS310128	3.0MM META SCREW - LONG THREAD 28MM
OS310130	3.0MM META SCREW - LONG THREAD 30MM
OS310132	3.0MM META SCREW - LONG THREAD 32MM
OS310134	3.0MM META SCREW - LONG THREAD 34MM
OS310612	3.0MM META SCREW - SHORT THREAD 12MM
OS310614	3.0MM META SCREW - SHORT THREAD 14MM
OS310616	3.0MM META SCREW - SHORT THREAD 16MM
OS310618	3.0MM META SCREW - SHORT THREAD 18MM
OS310620	3.0MM META SCREW - SHORT THREAD 20MM
OS310622	3.0MM META SCREW - SHORT THREAD 22MM
OS310624	3.0MM META SCREW - SHORT THREAD 24MM
OS310626	3.0MM META SCREW - SHORT THREAD 26MM
OS310628	3.0MM META SCREW - SHORT THREAD 28MM
OS310630	3.0MM META SCREW - SHORT THREAD 30MM
OS310632	3.0MM META SCREW - SHORT THREAD 32MM
OS310636	3.0MM META SCREW - SHORT THREAD 36MM
OS310640	3.0MM META SCREW - SHORT THREAD 40MM
OS310644	3.0MM META SCREW - SHORT THREAD 44MM

Twist off screws for Claw toe deformity & Metatarsalgia; Weil's Osteotomy	
Product Code	Description
OS310010	TWIST-OFF SCREW 10MM
OS310011	TWIST-OFF SCREW 11MM
OS310012	TWIST-OFF SCREW 12MM
OS310013	TWIST-OFF SCREW 13MM
OS310014	TWIST-OFF SCREW 14MM
OS310015	TWIST-OFF SCREW 15MM
OS310016	TWIST-OFF SCREW 16MM
OS310017	TWIST-OFF SCREW 17MM
OS310018	TWIST-OFF SCREW 18MM
OS309511	11MM TWIST-OFF SCREW 2.5MM CANCELLOUS THREAD
OS309512	12MM TWIST-OFF SCREW 2.5MM CANCELLOUS THREAD
OS309513	13MM TWIST-OFF SCREW 2.5MM CANCELLOUS THREAD
OS309514	14MM TWIST-OFF SCREW 2.5MM CANCELLOUS THREAD
OS309515	15MM TWIST-OFF SCREW 2.5MM CANCELLOUS THREAD
OS309517	17MM TWIST-OFF SCREW 2.5MM CANCELLOUS THREAD

Akin Osteotomy Staples	
Product Code	Description
OS310710	VARISATION STAPLE 90 DEGREE 10MM WIDE
OS310726	VARISATION STAPLE 26 DEGREE 10MM WIDE
OS310808	VARISATION STAPLE 90 DEGREE 8MM WIDE
OS310826	VARISATION STAPLE 26 DEGREE 8MM WIDE
OS310706	VARISATION COMPRESSION STAPLE 10MM WIDE

Foot Reconstruction System Disposables	
Product Code	Description
OS20870H	2.5MM META GUIDE WIRE 0.8MM X 70MM DOUBLE TROCAR
OS21070	K-WIRE DOUBLE TROCAR 1.0MM X 70MM
OS21070H	META GUIDE WIRE 1.0MM x 70MM DOUBLE TROCAR
OS292081H	2.5MM META GUIDE WIRE 0.8MM X 70MM TROCAR TIP
OS310101S	META SCREW CANN. DRILL LONG 2.0MM
OS310203	AO QC COMPRESSION SCREW DRILL 18MM
OS310204	AO QC COMPRESSION SCREW DRILL 12MM
OS310206	AO QC COMPRESSION SCREW DRILL 22MM
OS310304	2.5MM META SCREW DRILL 1.8MM
OS310305	2.5MM META SCREW COUNTERSINK DRILL



Ortho Solutions UK Limited  
West Station Business Park,  
Spital Road, Maldon, Essex.  
CM9 6FF, United Kingdom

Tel: +44 (0)1621 843 599  
Fax: +44 (0)1621 858 953  
Email: sales@orthosol.com

Ortho Solutions Inc - USA  
10901 West Toller Drive  
Suite #205, Littleton  
CO 80127

Tel: +1 (303) 495 5407  
Fax: +1 (303) 495 5408  
Email: us.sales@orthosol.com

Ortho Solutions Pty - Australia  
PO Box 1330,  
Milton,  
QLD. 4064, Australia

Tel: +61 (0)408 846821  
Email: salespty@orthosol.com

Rx only

Caution: Federal Law (USA)  
restricts this device to sale by  
or on the order of a physician

[www.orthosolutions.com](http://www.orthosolutions.com)

OS TD 00105\_18 - Rev 03 -  
Effective Date: Oct 2023

© 2023 Ortho Solutions UK Ltd. Ortho Solutions and the OS logo  
are trade marks of Ortho Solutions UK Ltd.



Consult  
instructions  
for use



<sup>2</sup> CE marking and the notified body number is applied per part number and  
appears on the device packaging, or the device if applicable.



Surgeon must be  
fully trained in the  
surgical technique

OrthoSolutions   
Group

Advancing Foot & Ankle Care