



Foot Osteotomy System

Product Rationale

OrthoSolutions 
Group

Advancing Foot & Ankle Care

Foot Osteotomy System Implants

Key features

The Foot Osteotomy System (FOS) can be used for various forefoot procedures including, but not limited to:

- Fixation of scarf and chevron osteotomies of the first metatarsal using the 3.0mm cannulated Scarf Screw
- Fixation of Akin or other phalangeal osteotomies using the 2.3mm cannulated Phalangeal Osteotomy Screw
- Fixation of Weil osteotomies and other lesser toe surgery procedures using the solid 2.0mm Weil Screw

Compact and simple instrumentation allows easy screw insertion

Headless, differential-pitch screws provide compression and reduce the incidence of post-operative soft-tissue irritation

OS Product code	Biomet™ equivalent code	Description
3.0mm Scarf Screw		
OS320010	V3010	Scarf Screw 10mm
OS320012	V3012	Scarf Screw 12mm
OS320014	V3014	Scarf Screw 14mm
OS320016	V3016	Scarf Screw 16mm
OS320018	V3018	Scarf Screw 18mm
OS320020	V3020	Scarf Screw 20mm
OS320022	V3022	Scarf Screw 22mm
OS320024	V3024	Scarf Screw 24mm
OS320026	V3026	Scarf Screw 26mm
OS320028	V3028	Scarf Screw 28mm
OS320030	V3030	Scarf Screw 30mm
OS320032	V3032	Scarf Screw 32mm
OS320034	V3034	Scarf Screw 34mm
OS320036	V3036	Scarf Screw 36mm
OS320038	V3038	Scarf Screw 38mm
OS320040	V3040	Scarf Screw 40mm



2.3mm Phalangeal Osteotomy Screw		
OS319008	Not Applicable	2.3mm Phalangeal Osteotomy Screw 8mm
OS319010	Not Applicable	2.3mm Phalangeal Osteotomy Screw 10mm
OS319012	V2312	2.3mm Phalangeal Osteotomy Screw 12mm
OS319014	V2314	2.3mm Phalangeal Osteotomy Screw 14mm
OS319016	V2316	2.3mm Phalangeal Osteotomy Screw 16mm
OS319018	V2318	2.3mm Phalangeal Osteotomy Screw 18mm
OS319020	V2320	2.3mm Phalangeal Osteotomy Screw 20mm
OS319022	V2322	2.3mm Phalangeal Osteotomy Screw 22mm
OS319024	V2324	2.3mm Phalangeal Osteotomy Screw 24mm
OS319026	V2326	2.3mm Phalangeal Osteotomy Screw 26mm
OS319028	V2328	2.3mm Phalangeal Osteotomy Screw 28mm




2.0mm Weil Screw		
OS320111	V2011	Weil Screw 11mm
OS320112	V2012	Weil Screw 12mm
OS320113	V2013	Weil Screw 13mm
OS320114	V2014	Weil Screw 14mm
OS320115	Not Applicable	Weil Screw 15mm



Foot Osteotomy System Instruments

Below illustrates kit listing for loan kits:



Number	Product code	Description	Quantity
1	OS320059	Sterilisation/instrument tray	1
2	OS320058	FOS bone holding forceps	1
3	OS320054	FOS depth gauge	1
4	OS320056	FOS drill guide	1
5	OS320150	FOS bone lever	1
6	OS320050	FOS AO QC handle	1
7	OS320053	Scarf screwdriver shank	1
8	OS319054	2.3mm phalangeal screw cannulated driver shank	1
9	OS319056	Weil screw solid driver shank	1

Sterile Consumables			
A*	OS319050	2.3mm phalangeal screw countersink	1
B*	OS320051	Scarf screw countersink	1
C*	OS320052	Scarf screw cannulated drill	1
D*	OS319052	2.3mm phalangeal screw cannulated drill	1
Not displayed	OS320049*	FOS solid 2.0mm drill	

*Only available sterile packed, products can be ordered in quantities of 10



Advancing Foot & Ankle Care



OrthoSolutions UK Limited,
West Station Business Park,
Spital Road, Maldon, Essex.
CM9 6FF, United Kingdom

Tel: +44 (0)1621 843 599
Fax: +44 (0)1621 858 953
Email: sales@orthosol.com

Ortho Solutions Inc - USA
10901 West Toller Drive
Suite #205, Littleton
CO 80127

Tel: +1 (303) 495 5407
Fax: +1 (303) 495 5408
Email: us.sales@orthosol.com

Ortho Solutions Pty - Australia
PO Box 1330,
Milton,
QLD. 4064, Australia

Tel: +61 (0)408 846821
Email: salespty@orthosol.com

www.orthosolutions.com

OS TD 00107_18 Rev 02
Effective Date: Sept 2024

© 2023 Ortho Solutions UK Ltd. Ortho Solutions and the OS logo
are trade marks of Ortho Solutions UK Ltd.



Consult
instructions
for use



2797²



Surgeon must be
fully trained in the
surgical technique

² CE marking and the notified body number is applied per part number
and appears on the device packaging, or the device if applicable.